

BIO+ HFC-02

Portable WIFI Hand-held Fundus Camera



www.portable-ultrasound-scanner

Quality, Professional, Innovational

Ergonomic Design

- 3.97 inch HD touch screen, 800*480
- 16 million pixels HD image
- Lightweight 560g,portable to take
- Feel comfortable and easy to operate
- Manual/Real time focus system

Android Intelligent System Design

- One hand operation with Android system
- DICOM interface
- 4G mobile network
- Work independently without any device

Wireless Transfer Design


- Support data storage,reading,wireless transfer,wireless print
- Real time storage through WIFI,BT,DICOM3.0,cloud,We-chat,Facebook, etc
- Wireless print via WIFI,BT and 4G

Product Parameter

Product category	Non-mydratic fundus camera
Product Screen	3.97 inch Multi-touch Screen, 800*480
Field of view	45°
Image Resolution	16MP
Target fixation	8 point internal targets fixation
Minimum pupil diameter	2.5mm
Flash illumination	White light&infrared
Focus	Manual/Auto
Focus range	-20D to + 20D
Image format	JEPG
Video format	3GP
Image color	Color,black and white
Memory	Internal 16 GB
Operation System	PC:Windows system APP:Android
Storage	USB2.0/WIFI/Bluetooth/4G
Connect function	USB2.0/WIFI/Bluetooth/4G/wireless printing etc
Image sharing	DICOM3.0,Cloud,Wechat,facebook,QQ,Microblog etc
Battery	3.7V rechargeable lithium battery
Adapter	50/60Hz,100~240V AC input;5V 1A DC output



Real Images



Wuxi Biomedical Technology Co., Ltd.

www.portable-ultrasound-scanner.com

Add: Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel: 86-510-85572789

Fax: 86-510-85572789

Email: sales@cnbiomedical.com

Mobile phone: 86-13382885910



is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.